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RETICULOCYTE STAINING AND COUNTING KIT (GB)

Bright cresyl and methylene blue stains precipitate with RNA (ribonucleic acid), the latter appearing as blue granules sometimes linked by a reticule.

Since the RNA of the reticulocytes disappears few days after they get into the blood, their counting will serve to know the number of cells released by bone marrow into blood.

The number of reticulocytes is probably the best and easiest indicator of erythropoiesis. An increase can be interpreted as a demand for new erythrocytes and proper functioning of the marrow. On the contrary, if the demand decreases or if the marrow does not function properly, then the number of reticulocytes decreases.

INSTRUCTIONS FOR USE

Add 3 drops of blood to a tube containing 100 µl of stable stain, mix well and incubate for 10-15 minutes at room temperature. Mix again, make the smear and let dry.

View through and oil immersion slide without further stain. In reticulocytes, the precipitated RNA is reticulate or granular and turns into intense blue. Erythrocytes turn into pale blue or greenish blue. Care must be taken since stain particles that adhere to the surface of the small corpuscles can be misinterpreted.

The coloration of preparations remains stable if they are not exposed to light.

The percentage of erythrocytes containing reticulate or granules is calculated with at least 1.000 erythrocytes in different parts of the extension must be observed. After calculating the percentage of reticulocytes, the total count of them can be determined multiplying the percentage for the erythrocytes count.

Normal results: adults usually have a reticulocyte count ranging between 0.5 and 2%. New-borns have a count ranging between 2.5 and 6.5%, decreasing to the adult range after 2 weeks.

RECOMMENDATIONS

Not suitable for any other application than its intended use. Do not use if packaging is damaged. Do not use if dirt is observed. Do not use the product if it has changed its original color. Do not use if you have exceeded the expiry date. Dispose as biological waste. Applicable safety data sheets can be requested. Storage: keep the product in a closed container, at room temperature, in a well-ventilated area, and away from any hot or ignition point.

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