

ACID PHOSPHATASE LEUKOCYTE

Cytochemical staining on blood or bone marrow smears in the classification of lymphoproliferative syndromes

10 x 4 tests

REF 3093

PREFACE

The kit has been designed to reduce the reagents volume and minimize the exposure of the operator to the chemicals, to simplify the procedure and the disposal of the reagents. Based on current knowledge, the least toxic and polluting reagents were used in the kit.

PRINCIPLE

Blood or bone marrow smears are incubated with Naphtol AS-BI phosphate and pararosaniline. A bright red color insoluble compound in cytoplasm indicates the presence of acid phosphatase in lymphocytes and monocytes. Staining also occurs in presence of tartrate. Intensity and frequency of stained granules in leucocytes and tartrate inhibition are evaluated under the microscope.

REAGENTS AND MATERIALS

Kit components:

*REAGENT 1 Sodium nitrite (lyophilized)	REF 3093
TOXICITY: Toxic substance. Do not swallow.	10 vials
*REAGENT 2 Pararosaniline	1 x 10 ml
TOXICITY: toxic substance if in contact with skin and inhalation. Keep away from light.	
REAGENT 3 Buffer	1 x 30 ml
REAGENT 4 Naphtol AS-BI phosphate	1 x 10 ml
PLATES Disposable multi well (4 wells in each plate)	10
COVER In black color for the plates	1

(* Dangerous reagents are market by an asterisk. Refer to MSDS.

STABILITY: sealed and stored at 2-8°C, reagents are stable up to the expiration date on the label.

ADDITIONAL REAGENT NOT PROVIDED WITH THE KIT

SODIUM TARTRATE

To classify hairy cells in which acid phosphatase is inhibited by tartrate. Determination of acid phosphatase with or without tartrate helps to evaluate the presence of an isoenzyme inhibited by such a compound in the cells.

REAGENTS REQUIRED BUT NOT PROVIDED

FIXATIVE:		
preparation:	formaldehyde 37%	1 volume
	absolute ethanol	9 volumes
COUNTERSTAINING:	Giemsa solution.	

MATERIALS REQUIRED BUT NOT PROVIDED

400x or 1000x **microscope** for slide reading.

Pipettes with 1 ml and 2.5 ml disposable tips or 1 ml and 3 ml graduated Pasteur pipettes for sampling and dispensing the reagents.

Thermostat set at 37°C.

Timer.

Deionized water.

SAMPLE

Blood (preferably from capillary) or bone marrow smears. Blood samples may be collected in EDTA or heparin.

Samples can be stored at room temperature (18-26°C) and protected from dust, for several days without any significant variation in activity.

Fixed slides can be stored for many weeks.

MANUAL ASSAY PROCEDURE

A) FIXATION OF THE SLIDES (see notes)

1. Fix the air-dried slides for 1 minute in the fixing solution.
2. Wash both sides of the slide in plenty of deionized water, drain it and wait till it is dry. The recommended fixative contains formaldehyde. Even a small quantity of aldehyde on the slides may inhibit the enzyme. Therefore, it is necessary to completely remove the fixative.

B) PREPARATION OF THE WORKING SOLUTION

Let the reagents reach room temperature before use. Unscrew the screw cap and carefully remove the rubber cap from a vial of Reagent 1.

1. Collect 1 mL of Reagent 2 using a pipette or Pasteur pipette and add it to a vial of Reagent 1. Replace the rubber cap and shake by inversion until the lyophilized reagent is completely dissolved. Wait for 2 minutes.

2. Reopen the vial and add 2.5 mL of Reagent 3.

3. Add 1 mL of Reagent 4, replace the cap, and shake well.

STABILITY: the working solution must be used right after it is prepared.

C) ACID PHOSPHATASE REACTION

1. Place the needed multi-well plates on a flat surface. Each plate and each bottle of working solution allows to run for 4 determinations.
2. Put the slides on the plate with the smear facing downwards, towards the bottom of the well, to ensure that the working solution come into contact with the smear.
3. Put the slide against one of the two long edges of the well. Between the longer side of the slide and the well, there will be a long groove into which the working solution will be injected.
4. Take 1 mL of working solution using a pipette or Pasteur pipette. Insert the tip into the central area of the opening and slowly inject the working solution. The solution will spread in the well, coming into contact with the smear. Less than 1 mL is enough to fill the well. Proceed the same way with the other slides.
5. Place the plate in a thermostat at 37°C and cover it with the black cover to protect it from light. If more plates are used, put them one over the other before covering. Incubate for 60 minutes.
6. Remove the slides with tweezers or fingers (wearing disposable gloves) and rinse them in running tap water. To facilitate this step, gently press one end of the slide so that the other one lifts up.
7. Washed and dried plates can be used for slide storage.

D) COUNTERSTAINING (see notes 1 and 2)

1. Counterstain with Giemsa solution for 10 minutes.
2. Rinse in running tap water, dry and read under the microscope. Experience in cytochemical techniques allows for the evaluation of slides without counterstaining.

RESULTS AND PATHOLOGY

The appearance of bright red granules in cell cytoplasm proves enzymatic activity.

This reaction helps to classify lymphoproliferative syndromes. In the lymphoid line, the single block positivity was related to T-phenotype. The intensity of the reaction in plasma cells of multiple myeloma was related to the prognosis of the disease. In hairy cells leukemia the reaction is not tartrate inhibited.

NOTES

Plates can be used for fixing and counterstaining the smear. In this case, put the slides as described in section C) and inject the fixing solution or the dye into the groove instead of the working solution. Proceed as described in sections A) and D) for fixation, counterstaining, and washes.

WASTE DISPOSAL

Dispose of reagents and used materials according to the regulations of your country.

BIBLIOGRAPHY

Available upon request.

MANUFACTURER



Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY
phone +39 045 6700870 website <http://www.fardiat.com>
e-mail: order@fardiat.com e-mail: fardiat@fardiat.com

KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalogue number
	temperature limits
	use by
	caution
	read instructions for use

