LYSOZYME STANDARD

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LYSOZYME STANDARD

Aqueous Lysozyme Standard (500 µg/ml)
for Lysozyme Measuring
in Urine, Serum and Plasma

1 x 1 ml

REF 7500

Additional available kit:

Lysozyme

REF 7104 – 7106

Kit for turbidimetric determination of Lysozyme in urine, serum and plasma.

PRINCIPLE

The Lysozyme Standard solution contains 500 µg/ml of chicken egg lysozyme with activity equal to 50 KU/ mg of protein.

To draw the calibration curve for the lysozyme kit, dilute the Standard Lysozyme solution with saline solution to obtain working solutions with lysozyme concentrations between 2 and 50 µg/mL.

For example, dilute the Lysozyme Standard solution as it follows:

Standard		Saline solution		Lysozyme concentration
0.200 ml	+	1.80 ml	=	50 µg/ml
0.100 ml	+	1.15 ml	=	40 µg/ml
0.120 ml	+	1.88 ml	=	30 µg/ml
0.100 ml	+	2.40 ml	=	20 µg/ml

For concentrations lower than 20 µg/ml, dilute the Solution 50 µg/ml as it follows:

Soluzione 50 µg/ml		Saline solution		Lysozyme concentration
0.200 ml	+	0.8 ml	=	10 µg/ml
0.100 ml	+	0.9 ml	=	5 µg/ml
0.100 ml	+	2.4 ml	=	2 µg/ml

REAGENTS

The reagent is an aqueous standard ready to use.

Kit component:

Lysozyme standard (500 µg/ml)

REF 7500

1 x 1 ml

STABILITY: stored at 2-8°C, the reagent is stable up to the expiration date on the label.

MANUAL ASSAY PROCEDURE

Follow the determination procedure described in the lysozyme kit instruction sheet, using **8 µl** of **LYSOZYME STANDARD properly diluted** instead of the sample.

To draw the calibration line, note in a graphic the ΔAbs (2 min) obtained for the different Lysozyme Standard dilutions for the lysozyme concentration.

NOTES

- (*) Dangerous reagents are marked by an asterisk. Refer to MSDS.
- If the ΔAbs (2 min) obtained with the Lysozyme Standard solution diluted at 50 µg/ml is lower than 0.1 spectrophotometric unit, define the activity with higher lysozyme concentrations.
- If temperatures or wavelengths different than the ones stated are used, adjust properly also the concentration interval of the solutions.
- Disposal all waste material in accordance with the law.

REFERENCES

- Rudders et Bloch. Am.J.medical sciences (1971: 79-85)
- Horpcsy et al. Clin. Chem. 24/1 74-79. 1978)

IVD

CE

Issue 02 - Mar 2015 RR

Manufactured by: **FAR** srl

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

Phone +39 045 6700870 - Fax +39 045 7157763

website: <http://www.farddiag.com> e-mail: farddiag@farddiag.com