ACE CALIBRATOR

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ACE CALIBRATOR

Solutions for Wavelength Calibration and Spectrophotometer Sensitivity for Angiotensin Converting Enzyme Determination in Serum

1 x 2 ml

REF 7512

PRINCIPLE

For ACE determination, the spectrophotometric reading is performed in a spectrophotometer area where even a small wavelength variation defines a considerable change of the extinction coefficients of the substrate and of the product, and therefore of the $\Delta\epsilon$ employed to calculate the enzymatic activity. This kit allows the correction of this factor for the equipment used.

REAGENTS

Kit components:

REF 7512

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furanacryloylphenylalaninglycylglycine in buffered solution

1 vial

REAGENT P (Product)

furanacryloylphenylalanine in buffered solution

1 vial

STABILITY: stored at 2-8°C in its original package, the reagent is stable up to the expiration date on the label.

MANUAL ASSAY PROCEDURE

Wavelength: spectrophotometer 340 nm
Optical path: 1 cm
Reading: against distilled water
Temperature: 25°C

Read the REAGENT S (As) and REAGENT P (Ap) absorbencies.

Perform the following calculation:

$$\Delta\epsilon_c = (A_s - A_p) \times 2$$

CORRECTION OF THE CALCULATION FACTOR

Multiply the calculation factor used for ACE-S activity in U/L, amounting to 863, by 0.85 and divide by the $\Delta\epsilon_c$ value found. Therefore the calculation formula is the following:

$$\text{ACE activity in U/L} = (A_1 - A_2) \times 863 \times 0.85 / \Delta\epsilon_c$$

NOTES

- (*) Dangerous reagents are marked by an asterisk. Refer to MSDS.
- The solutions may also be employed to measure the $\Delta\epsilon$ value following J.E. Buttery and B.R. Chamberlain's method. Clin. Chem. 31 (1985) 645.

REFERENCE

- Harjanne A. Clin. chem. 30 (1984) 901

Issue 02 - Mar 2015 RR

Manufactured by: FAR srl

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY
Phone +39 045 6700870 - Fax +39 045 7157763
website: <http://www.fardiag.com> e-mail: fardiag@fardiag.com