

Direct BILIRUBIN

DPD Colorimetric method

R1: 2 x 100 ml + R2: 2 x 25 ml

CL09-250

INTENDED USE

Kit for quantitative determination of Direct Bilirubin in serum.

CLINICAL MEANING

Analysis of total and direct bilirubin is made to determine the presence of hepatic damages or diseases e.g. obstruction of bile ducts, haemolytic anemias, metabolic issues, stones. A typical sign of high bilirubin levels is jaundice, which manifests with yellow skin and yellow sclerae. An increase in direct bilirubin is a symptom of a hepatic problem, especially of the biliary excretion tract.

PRINCIPLE

In an acid medium, total bilirubin reacts with 3,5-DPD forming a red-colored diazo compound (azobilirubin), whose intensity is proportional to the concentration of direct bilirubin present in the sample.

While direct bilirubin, that is conjugated with glucuronic acid, is water-soluble and reacts directly, total bilirubin is obtained through the presence of an accelerator that splits the bond with albumin.

Indirect bilirubin can be calculated by the difference between total bilirubin and direct bilirubin.

SAMPLE

Non hemolyzed serum.

Analyze samples within 2 hours after collection. Protect samples from light.

STABILITY: 12 hours in the fridge at 2-8°C, 3 months at -20°C if protected from light.

REAGENTS

Only for in Vitro diagnostic use. Liquid mono-reagent ready to use.

Package content		CL09-250
REAGENT 1 Hydrochloric acid 170 mM		CL09-250R1 2 x 100 ml
REAGENT 2 3,5-DPD 0,5 mM Hydrochloric acid 170 mM		CL09-250R2 2 x 25 ml

Stability: Store at 2-8°C and protect from light to keep the reagents stable up to the expiration date on the label. Keep bottles closed when not in use. Do not use turbid reagents.

NECESSARY ITEMS – NOT PROVIDED

Usual laboratory equipment: UV/VIS Spectrophotometer with temperature control; automatic micropipettes; Optical glass cuvettes or, alternatively, disposable ones in optical polystyrene; Saline solution.

MANUAL ASSAY PROCEDURE

Primary wavelength:	546 nm
Reference wavelength:	800 nm
Optical Path:	1 cm
Reading:	Against blank sample
Temperature:	37 °C
Method:	Increasing endpoint
Reaction:	two final points
Ratio:	2/20/5

Bring reagents to the chosen temperature for the analysis.

Pipette in cuvette:

Pipette in cuvette:

	Calibrator	Sample
Reagent 1	1 ml	1 ml
Sample	---	100 µl
Calibrator	100 µl	---

Stir carefully. After exactly 5 minutes of incubation at 37°C, read the Blank Calibrator (AbsBCal) and the blank sample (AbsBS) against water.

Add:

Reagent 2	250 µl	250 µl
-----------	--------	--------

Stir carefully. After exactly 5 minutes of incubation at 37°C, read the Calibrator (AbsCal) and the Sample (AbsS) against water

Reaction volumes can be proportionally varied without any change in calculation.

CALCULATION

Calculate the concentration in the sample using the following formula:

AbsS – AbsBS

_____ x Val. Cal = BIL-D (mg/dl)

AbsCal – AbsBCal

REFERENCE VALUES

Up to 0.35 mg/dl (5.1 µmol/L).

Each laboratory should define its own reference values for this method.

QUALITY CONTROL – CALIBRATION

All Clinical Chemistry laboratories should implement a quality control program. Control serums of human origin are available for this purpose on request:

PRE-NORM serums with normal values

PRE-PATH serums with pathological values

If the method requires it, a multi-parameter calibrator of human origin is available (REF 7532).

PERFORMANCE CHARACTERISTICS

Sensitivity: the sensitivity of the method is 0.06 mg/dl.

Linearity: up to 20 mg/dl.

For higher values, dilute the samples with saline solution and multiply the result by the dilution factor.

Precision:

Within run (n=30)	Average [mg/dl]	CV %
Sample 1	0.21	2.8
Sample 2	3.12	0.9

Between run (n=20)	Average [mg/dl]	CV %
Sample 1	0.24	3.12
Sample 2	4.8	1.4

Interferences:

- up to 500 mg/dl of hemoglobin does not interfere.

- up to 1000 mg/dl triglycerides do not interfere.

Direct light can cause a decrease of direct bilirubin up to 50% in an hour.

Correlation against a reference method: the correlation of FAR method (Y) against a reference method (X) gives this equation:

$$Y = 0.991x - 0.100 \quad (n=30) \quad r = 0.997 \quad (n=30)$$

DISPOSAL

The product must be used for professional analysis only. The product must be disposed of according to national/international laws.

WARNINGS AND PRECAUTIONS

The reagents may contain non-reactive components and various preservatives. Contact with the skin and ingestion should be avoided. Use the normal precautions expected with correct behaviour in laboratory.

REFERENCES

1. Pearlman F.C., Lee R.T.Y., Clin. Chem. 20, 447, (1974)
2. Blumenfeld T.A. et al., Am. J. Clin. Path. 69, 388 (1978)

MANUFACTURER

FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY








tel +39 045 6700870

website <http://www.farddiag.com>

e-mail: order@farddiag.com

e-mail: farddiag@farddiag.com

KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalog number
	temperature limits
	use by
	caution
	consult accompanying documents

Issue 01 - Jan 2021 RR