

CHOLESTEROL HDL

Direct Method - Selective Detergent

R1: 2 x 60 ml + R2: 2 x 20 ml

CL24-160

INTENDED USE

Kit for quantitative determination of Cholesterol HDL in serum and plasma.

CLINICAL MEANING

Cholesterol is a fat present in the blood, involved in various processes that are fundamental for the functioning of the body. If present in excessive quantities, it is one of the major risk factors for heart disease. Cholesterol is transported within molecular structures known as lipoproteins, more specifically low-density lipoproteins (LDL) and high-density lipoproteins (HDL).

HDL-C is mainly made by proteins and by a reduced quantity of cholesterol. It's considered as protective because it collects the excess cholesterol and carries it to the liver for removal. For this reason it is also called "good cholesterol".

PRINCIPLE

LDL, VLDL and Chylomicron (CM) react with PVS (polyvinyl sulfonic acid) and PEGME (polyethylene-glycol methyl ether) and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHOD) and cholesterol esterase (CHER). The enzymes selectively react with HDL to produce H₂O₂ which is quantified by the Trinder reaction.

SAMPLE

Serum, plasma (EDTA, Citrate, Li Heparine)

Use fresh samples.

REAGENTS

Only for in Vitro diagnostics. Liquid monoreagent ready to use.

| Package content | CL22-160 | Quantity |
|--|------------|-----------|
| REAGENT 1 MES Buffer, TODB, PVS, PEGME, detergent, EDTA. | CL22-160R1 | 2 x 60 ml |
| REAGENT 2 MES Buffer, Cholesterol Esterase, Cholesterol Oxidase, Peroxidase, 4-aminoantipyrine, detergent. | CL22-160R2 | 2 x 20 ml |

Stability: reagents are ready to use. Store at 2-8°C and protect from light to keep the reagents stable up to the expiration date on the label. Do not freeze. Once opened the reagents are stable for 30 days at 2-8°C if contamination is avoided. Keep bottles closed when not in use. Do not use turbid reagents.

NECESSARY ITEMS – NOT PROVIDED

Usual laboratory equipment: UV/VIS Spectrophotometer with temperature control; automatic micropipettes; Optical glass cuvettes or, alternatively, disposable ones in optical polystyrene; Saline solution.

LDL Calibrator

MANUAL ASSAY PROCEDURE

| | |
|-----------------------------|-----------------------|
| Analysis: | increasing endpoint |
| Wavelength: | 600 nm |
| Cuvette: | 1 cm optical path |
| Temperature: | 37°C |
| Reaction Time: | 5 + 5 minutes |
| Reading: | against blank reagent |
| Sample/Reagent 1/Reagent 2: | 1/75/25 |

Let reagents necessary to perform the test reach the chosen temperature for the analysis.

Pipette in cuvette:

| | Blank Reagent | Calibrator | Sample |
|-----------------|---------------|------------|--------|
| Distilled water | 3 µl | - | - |
| Calibrator | - | 3 µl | - |
| Sample | - | - | 3 µl |
| Reagent 1 | 225 µl | 225 µl | 225 µl |

Mix and incubate for 5 minutes at 37°C.

Read the absorbance of the calibrator (AbsC 1) and the sample (AbsS 1) against the blank reagent. Then pipette:

| | | | |
|-----------|-------|-------|-------|
| Reagent 2 | 75 µl | 75 µl | 75 µl |
|-----------|-------|-------|-------|

Mix and incubate for 5 minutes at 37°C. Then read the absorbance of the calibrator (AbsC2) and the sample (AbsS 2) against the blank reagent.

Reaction volumes can be proportionally varied without any change in calculation.

CALCULATION

Calculate the concentration in the sample using the following formula:

[mg/dl] HDL- cholesterol =

$$(\text{AbsS2} - \text{AbsS1}) / (\text{AbsC2} - \text{AbsC1}) \times \text{calibrator conc.}$$

[mmol/L] HDL-cholesterol =

$$(\text{AbsS2} - \text{AbsS1}) / (\text{AbsC2} - \text{AbsC1}) \times \text{calibrator conc.} \times 0.02586$$

REFERENCE VALUES

The expected values for HDL Cholesterol are as follow:

< 40 mg/dl (1,03 mmol/L)

40– 59 mg/dl (1,03-1,53 mmol/L)

≥ 60 mg/dl (≥1.55 mmol/L)

a major risk factor for heart diseases

the higher uour HDL, the better

protective against heart disease

Each laboratory should define its own reference values for this method.

QUALITY CONTROL – CALIBRATION

A quality control program is recommended for all clinical laboratories.

Control serums of human origin are available for this purpose on request:

PRE-NORM serums with normal values

PRE-PATH serums with pathological values

If the method requires it, a multi-parameter calibrator of human origin is available.

PERFORMANCE CHARACTERISTICS

Sensitivity

The sensitivity of the method is 1,1 mg/dl.

Linearity: up to 180 mg/dl.

For higher values, dilute the samples 1:5 with saline solution and multiply the result by 5.

Precision

| Within run (n=20) | Average [mg/dl] | SD | CV % |
|-------------------|-----------------|------|------|
| Sample 1 | 29,0 | 0,3 | 1,0 |
| Sample 2 | 53,07 | 0,41 | 0,8 |
| Sample 3 | 90,56 | 0,84 | 0,9 |

| Between-run (n=20) | Average [mg/dl] | DS | CV % |
|--------------------|-----------------|------|------|
| Sample 1 | 29 | 0,65 | 2,3 |
| Sample 2 | 53,07 | 1,36 | 2,02 |
| Sample 3 | 90,56 | 2,02 | 2,2 |

Interferences: up to 40 mg/dl of bilirubin does not interfere. Up to 1000 mg/dl of hemoglobin does not interfere. Up to 10 mM of ascorbic acid does not interfere. Up to 1000 mg/dl of lipids (as triglycerides) does not interfere.

Correlation against a reference method: the correlation of the method (Y) against a reference method (X) gives this equation:

$$Y = 1,05X - 4,69$$

$$r = 0,986$$

DISPOSAL

The product must be used for professional analysis only. The product must be disposed of according to national/international laws.

WARNINGS AND PRECAUTIONS

The reagents may contain non-reactive components and various preservatives. Contact with the skin and ingestion should be avoided. Use the normal precautions expected with correct behaviour in laboratory.

REFERENCES

- Castelli WP e al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
- PisaniT, Gebski CP, Leary Et, et al. Accurate Direct Determination od Low-Density Lipoprotein Cholesterol Assay. Arch Pathol KLab Med 1995; 119:1127.

MANUFACTURER

FAR








Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

tel +39 045 6700870

website <http://www.farddiag.com>

e-mail: order@farddiag.com e-mail: farddiag@farddiag.com

KEY SYMBOLS

| | |
|---|------------------------------------|
|  | In Vitro diagnostic medical device |
|  | batch number |
|  | catalogue number |
|  | temperature limits |
|  | use by |
|  | caution |
|  | consult accompanying documents |

Issue 01 - Jan 2021 RR