

# PHENOLS

Colorimetric quantitative determination of phenolic compounds on urine

50 tests

**REF** CM07-50T

## INTENDED USE

Kit for quantitative *in vitro* determination of Phenolic Compounds on urine.

## PRINCIPLE

In alkaline medium, phenolic compounds react with a 4-aminoantipyrine condense. The developed product is oxidized with ferricyanide, forming a red color complex which can be measured photometrically.

## PURPOSE OF THE TEST

The determination of urinary phenolic compounds is a useful diagnostic assay for patients with supposed hyperthyroidism, mellitus diabetes caused by hyperthyroidism, tumors producing catecholamines and its metabolites. This method can also be used as a screening test in supposed benzene intoxication.

## REAGENTS

Kit components:

<b>REAGENT 1</b> 4-aminoantipyrine	<b>REF</b> CM07-50T CM07-50R1: 1 x 80 ml
<b>REAGENT 2</b> Ferricyanide	CM07-50R2: 1 x 80 ml
* <b>REAGENT 3</b> Extraction solution	CM07-50R3: 2 x 90 ml
<b>STANDARD</b> Phenol 500 mg/L	CM07-50TS: 1 x 3 ml

(\* Dangerous reagents are marked by an asterisk. Refer to MSDS.

**STABILITY:** stored at 2-8°C, sealed reagents are stable up to the expiration date on the label.

## REQUIRED EQUIPMENT

Centrifuge, spectrophotometer or filter photometer at 460 nm (450 - 500 nm).

## DILUTION OF THE STANDARD

Dilute the standard 1:10 with distilled water, to get a standard with a 50 mg/L concentration.

**STABILITY:** 1 month at 2-8°C.

## SAMPLE

24-hour urine.

**NOTICE:** store the urine at 2-8°C and perform the assay as soon as possible.

## MANUAL ASSAY PROCEDURE

Wavelength:	460 nm (450 – 500 nm)
Optical path:	1 cm
Reading:	against blank reagent
Temperature:	room temperature
Method:	colorimetric endpoint
Linearity:	up to 200 mg/L
Minimum sensitivity:	0.1 mg/L
C.V. (intra-assay):	2.77 %

Pipette into 3 centrifuge tubes labeled as follows:

	Blank Reagent	Sample	Standard
Distilled water	0.2 ml	---	---
Sample	---	0.2 ml	---
Diluted standard	---	---	0.2 ml
Reagent 1	1.5 ml	1.5 ml	1.5 ml

Mix thoroughly and add:

Reagent 2	1.5 ml	1.5 ml	1.5 ml
-----------	--------	--------	--------

Mix carefully and let it stand for 3 minutes. Then add:

Reagent 3	3.5 ml	3.5 ml	3.5 ml
-----------	--------	--------	--------

Shake and centrifuge at 2500-3000 rpm for 5 minutes. Read the supernatant absorbances (Ac e Astd) against the blank reagent.

## CALCULATION

Phenols (mg/L) = (Ac / Ast) x 50

Phenols mg/litro x L of 24 h urine = Phenols mg/24 hours

## REFERENCE VAUES

Up to 50 mg/24 hours.

## PERFORMANCE CHARACTERISTICS

**Sensitivity:** the sensitivity of the method is 0.1 mg/L.

**Linearity:** up to 200 mg/L.

For higher values, dilute the sample 1:10 with saline solution and multiply the result by 10.

**Precision:**

Between run (n=10)	Mean [U/L]	CV %
Sample 1	300	2.9
Sample 2	100	2.7

**Correlation against a reference method:** the correlation of FAR method (Y) against a reference method (X) gives a correlation of 0.9856

## DISPOSAL

The product must be used for professional assay only. Dispose of the product according to national/international laws.

## WARNINGS AND PRECAUTIONS

**REAGENT 1 and REAGENT 2:** Not dangerous



**STANDARD**

Advisory: Danger

Hazard Statements:

**H315** Causes skin irritation.

**H319** Causes severe eyes irritation.

**H335** May cause respiratory irritation.

## NOTE

Comparing FAR kit (y) with Yamaguci and Hayashi method to define phenols, the obtained correlation coefficient is 0.985.

## REFERENCE

1. Y. Yamaguchi et C. Hayashi, Clin. Chem., 23/11, 2151-2154 (1977)

## KEY SYMBOLS

<b>IVD</b>	In Vitro diagnostic medical device
<b>LOT</b>	batch number
<b>REF</b>	catalogue number
	temperature limits
	use by
	caution
	read instructions for use

**IVD**

**CE**

Ed. 01 - July 2022

## MANUFACTURER



FAR  
Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY  
phone +39 045 6700870  
website <http://www.fardiag.com>  
e-mail: [order@fardiag.com](mailto:order@fardiag.com) e-mail: [fardiag@fardiag.com](mailto:fardiag@fardiag.com)