

Total LIPIDS liquid

Colorimetric determination
on serum

4 x 100 ml

REF CM09-400

PRINCIPLE

Total lipids in serum, after a heat treatment with sulphuric acid, react with phosphoric acid-vanillin mixture, developing a pink color whose intensity is proportional to the lipids concentration.

REAGENTS

Kit components:

(*) REAGENT 1 REF CM09-400 Quantity CM09-400R1 4x100 ml

Vanillin 7.8 mmol/L

Phosphoric acid 12 mmol/L

STANDARD (Std) CM09-400S 1x4 ml

Cholesterol 800 mg/dl (26.68 mmol/L)

STABILITY: stored at 15-30°C and tightly closed the reagents are stable up to the expiration date on the label

ADDITIONAL REAGENT NOT INCLUDED WITH THE KIT

Ultra pure sulphuric acid d = 1.84

SAMPLE

Serum

MANUAL ASSAY PROCEDURE

Wavelength: 530 nm (510-550 nm)

Optical path: 1 cm

Reading: against blank reagent

Temperature: 100°C for lipo-sulphuric mixture
room for the reaction

Method: endpoint

Reaction: 20 minutes

Linearity: up to 1500 mg/dl

Sample/reagent: 1/25

Pipette into pipex tubes labeled as follows:

	Sample	Standard
Sulphuric acid (d = 1.84)	2500 µl	2500 µl
Standard	---	100 µl
Sample	100 µl	---

Mix accurately and incubate in a hot bain-marie for 10 minutes.

Cool.

These are the lipo-sulphuric reaction mixtures of both sample and standard to be used to run the test.

Pipette into tubes or cuvettes labeled as follows:

B/R: blank reagent, S: sample, Std: standard:

	B/R	S	Std
Reagent 1	2500 µl	2500 µl	2500 µl
Standard lipo-sulphur. mixture	---	---	100 µl
Sample lipo-sulph. mixture	---	100 µl	---
Sulphuric acid (d = 1.84)	100 µl	---	---

Mix the tubes well and incubate for 20 minutes at 20-25°C.

Read the standard (Astd) and sample (As) absorbances against the blank reagent, within 45 minutes.

CALCULATION

lipids (mg/dl) = (As / Astd) x 800

lipids (mmol/L) = (As / Astd) x 20.68

REFERENCE VALUES

Lipids in serum: 450-800 mg/dl

QUALITY CONTROL – CALIBRATION

A quality control program is recommended for all clinical chemistry laboratories.

PERFORMANCES CHARACTERISTICS

Linearity: up to 1500 mg/dl.

For higher values, dilute the samples and multiply the result by the dilution factor.

Within-run precision:

	Level 1	Level 2
Mean (mg/dl)	532	1105
DS	3.90	9.32
CV %	0.73	0.84

Between-run precision:

	Level 1	Level 2
Mean (mg/dl)	552	980
DS	6.00	19.5
CV %	1.09	1.99

Correlation

FAR kit for total lipids determination shows a correlation coefficient of 0.987 compared to another kit available on the market.

DISPOSAL

The product must be used for professional assay only. Dispose of the product according to national/international laws.

WARNINGS AND PRECAUTIONS



REAGENT 1

H314 Causes severe skin burns and eyes damage.



STANDARD

H225 Highly flammable liquid and vapor.

NOTES

- (*) Dangerous reagents are marked by an asterisk. Refer to MSDS.
- Reaction volumes can be proportionally changed.
- Use only glassware without any trace of detergents.
- Reagents require an accurate pipette due to their viscosity.
- Chemistry analyzer parameters are available.

REFERENCE

- Zöllner N., Kirsch K., Z. ges. exp. Med., 135, 545 (1962)

MANUFACTURER



FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

tel +39 045 6700870

website <http://www.farddiag.com>

e-mail: order@farddiag.com

e-mail: farddiag@farddiag.com



KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalog number
	temperature limits
	use by
	caution
	read instructions for use

Ed. 01 - Jan 2022