

ACE

UV Determination of Angiotensin Converting Enzyme in Serum and Plasma

9 x 4 mL**REF CY02-36**

Other available kits:

ACE (liquid reagents)**REF CY02-36**

Available for quality control:

ACE-CONTROL SERUM N + P**REF 7508**

Control serums in normal and pathological range

ACE-CALIBRATOR**REF 7512**

For an accurate control of instrument calibration

ACE-STANDARD**REF 7511**

Standard of ACE for Measuring of the Enzyme in Serum

PRINCIPLE

The angiotensin converting enzyme (ACE) catalyzes hydrolysis of furylacryloylphenylalanyl-glycylglycine (FA-Phe-Ala-Gly-Gly) substrate to furylacryloyl phenyl-alanine and glycylglycine.

Hydrolysis is related to an absorbance decrease valued at 340 nm and is proportional to enzymatic activity.

REAGENTS

Kit composition:

REF CY02-36 **Quantity****REAGENT 1** (lyo)
FA-Phe-Ala-Gly-Gly**7004R1** **9 vials****REAGENT 2**
Buffer pH 8,4**7004R2** **2 x 20 mL**

STABILITY: stored at 2-8°C, reagents are stable up to expiration date.

PREPARATION OF WORKING REAGENT

Reconstitute the contents of a vial of Reagent 1 with exactly 4.2 mL of Reagent 2. Shake gently until complete dissolution.

STABILITY: 20 days at 2-8°C away from light.

SAMPLE

Serum or heparinized plasma.

STABILITY: 4 days at , 6 months at -20°C.

MANUAL ASSAY PROCEDURE

Wavelength: 340 nm
Optical path:: 1 cm
Reading: against distilled water
Temperature: 37°C
Method: fixed time
Reaction: 15 minutes
Linearity: up to 250 U/L
Sample/Reagents: 1/10

NOTE: spectrophotometric reading is made in a substrate spectrum zone where even a small wavelength change corresponds to a high variation of extinction coefficients.

For a proper use, carefully check the wavelength calibration and the instrument sensibility.

For this aim, use the product ACE CALIBRATOR.

Pipette into cuvette:

Working reagent	1.0 mL
Sample	0.1 mL

Mix and incubate at 37°C. After 5 minutes read A1 absorbance and after exactly 15 minutes from the first reading, read A2 absorbance.

CALCULATION

ACE Activity (in U/L): = (A1-A2) x 863

Reaction volumes can be proportionally varied without any change in calculation. Each laboratory should define its own reference values for this method.

REFERENCE VALUES

AVERAGE	±	DS
90.1 U/L	±	24.3 U/L

QUALITY CONTROL – CALIBRATION

All Clinical Chemistry laboratories should implement a quality control program.

Control serums of human origin are available for this purpose on request:

ACE CONTROL SERUM N+P with normal and pathological value ranges

A Standard for an accurate control is also available:

ACE STANDARD 2x1 ml**PERFORMANCE CHARACTERISTICS**Linearity: up to 250 u/l

For higher values, dilute the sample 1:2 with saline solution, repeat the test and multiply the result by 2.

Within-precision:

	Level 1	Level 2
Average [[U/L]	75.3	153
DS	0.74	2.85
CV %	0.98	1.86

Between-precision:

	Level 1	Level 2
Average [U/L]	81.0	150
DS	1.45	4.59
CV %	1.79	3.06

Correlation:

FAR ACE kit shows a correlation coefficient equal to 0.987 in comparison to another kit available on the market.

NOTES

- Refer to MSDS.
- ACE is a metal-protein, hence avoid chelates in the sample preparation.
- Reaction volumes can be proportionally varied.
- To calculate activity use the following formula:

$$U/L = (A1 - A2) \times [(Vt \times 1000) / (\Delta\epsilon \times l \times Vs \times t)]$$
 where:
 A1: absorbance in the sample after 5 minute incubation;
 A2: absorbance in the sample after 15 minute incubation from the first reading;
 Vt: total volume (reagent +sample) in ml;
 $\Delta\epsilon$: variation of extinction coefficient at 340 nm;
 l: optical path in cm;
 Vs: sample volume in ml;
 t: incubation time in minutes.

Under these test conditions the formula becomes:

$$U/L = (A1 - A2) \times [(1,1 \times 1000) / (0,85 \times 1 \times 0,1 \times 15)] =$$

$$= (A1 - A2) \times 863$$

$\Delta\epsilon$ was defined by research spectrophotometers. Using chemical analyzers, $\Delta\epsilon$ might reach a different value with following modification of U/L values in normal and pathologic people. Use ACE CALIBRATOR to calculate $\Delta\epsilon$ for the instrument used.

- Each laboratory should define its own reference values.
- Chemistry analyzer parameters are available.

DISPOSAL

The product must be used for professional analysis only. The product must be disposed of according to national/international laws.

WARNINGS AND PRECAUTIONS

The reagents may contain non-reactive components and various preservatives. Contact with the skin and ingestion should be avoided. Use the normal precautions expected with correct behaviour in laboratory.

REFERENCES

- Harjanne A. Clin. Chem. 30 (1984) 901

MANUFACTURER

FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

tel +39 045 6700870

website <http://www.fardiag.com>e-mail: order@fardiag.come-mail: fardiag@fardiag.com**KEY SYMBOLS**