


Kinetic Determination of enzyme
Glucose 6-Phosphate Dehydrogenase
in Serum, Plasma and Erythrocytes

4 x 25 ml

REF CY06-100

PRINCIPLE

The activity of the enzyme Glucose-6-Phosphate Dehydrogenase (E.C. 1.1.1.49) is defined by measuring the speed of the absorbance increase at 340 nm due to the reduction of NADP⁺ according to the following reaction:



REAGENTS

Kit components:

REAGENT 1

Triethanolamine buffer pH 7.6
EDTA, sodium azide 15 mmol/L

REAGENT 2 (lyo)

NADP⁺
Glucose-6-P

*REAGENT 4

Digitonine

STABILITY: store at 2-8°C to keep the reagents stable up to the expiration date on the label.

REF CY06-100
CY06-100R1 **Quantity**
4 x 25 ml

CY06-100R2 **1 vial**

CY06-100R3 **1 vial**

CY06-100R4 **1 x 20 ml**

PREPARATION OF WORKING REAGENTS

REAGENT 2

Dissolve the contents of the vial of Reagent 2 with exactly 3.0 ml of distilled water. Shake gently until complete dissolution.

STABILITY: 4 weeks at 2-8°C.

REAGENT 3

Dissolve the contents of the vial of Reagent 3 with exactly 1.5 ml of distilled water. Shake gently until complete dissolution.

STABILITY: 4 weeks at 2-8°C.

SAMPLE

Serum or plasma: use fresh and absolutely hemolysis free samples.
Erythrocytes (hemolized): follow Procedure 2 to prepare the hemolized.

MANUAL ASSAY PROCEDURE

Wavelength: Hg 365 nm, 340 nm or Hg 334 nm
Optical path: 1 cm
Reading: against air
Temperature: 25°C, 30°C, 37°C (serum or plasma)
25°C (erythrocytes)
Reaction time: 3 minutes

1. MEASURE OF G6P-DH IN SERUM OR PLASMA

Let the reagents reach the chosen temperature for the analysis.

Pipette in cuvette:

	25°C – 30°C	37°C
Reagent 1	1.00 ml	1.00 ml
Reagente 2	50 µl	50 µl
Serum/plasma	500 µl	250 µl

Mix, incubate for 10 minutes at the test temperature. Add:

Solution 3	25 µl	25 µl

Mix. Read initial absorbance and start timer simultaneously. Repeat the reading after exactly 1, 2 and 3 minutes. Determine the average value of absorbance change per minute (ΔA/min).

CALCULATION

Calculate the G6P-DH activity (mU/ml) in serum or plasma using the following formulas:

	340 nm	Hg 365 nm	Hg 334 nm
25°C/30°C	500 x ΔA/min	900 x ΔA/min	510 x ΔA/min
37°C	841 x ΔA/min	1514 x ΔA/min	858 x ΔA/min

2. MEASURE OF G6P-DH IN ERYTHROCYTES

Define the number of erythrocytes per ml of blood as normal procedure.

PREPARATION OF SAMPLE

Wash 0.2 ml of blood for three times with 2 ml of saline solution. After each washing, centrifuge for 10 minutes at approximately 3000 rpm/minute.

Re-suspend the washed erythrocytes and centrifuged in 0.5 ml of Reagent 4, incubate for 15 minutes at approximately 4°C, then centrifuge again. To define the activity, use the supernatant. Perform the test within 2 hours of the preparation of the hemolized.

Let the reagents reach the chosen temperature for the analysis.

Pipette in cuvette:

Reagent 1	1.00 ml
Reagent 2	30 µl
Hemolized	15 µl

Mix, incubate for 5 minutes at 25°C. Add:

Reagent 3	15 µl
-----------	-------

Mix. Read initial absorbance and start timer simultaneously. Repeat the reading after exactly 1, 2 and 3 minutes. Determine the average value of absorbance change per minute (ΔA/min).

CALCULATION

To calculate the G6P-DH activity in erythrocytes use the following formulas:

	340 nm	Hg 365 nm	Hg 334 nm
mU (erythrocytes in 1 ml of blood*)	33650 x ΔA/min	60571 x ΔA/min	34304 x ΔA/min

* The G6P-DH activity is expressed as mU/10⁹ erythrocytes or as mU/g hemoglobin.

To calculate the G6P-DH activity as mU/10⁹ erythrocytes, as it is required for a comparison with the normal value, divide the calculated activity (mU in the erythrocytes present in 1 ml of blood) by the number of erythrocytes per ml of blood with erythrocytes count. For example, if you obtain 695 mU for the erythrocytes present in 1 ml of blood, and the result of the erythrocytes count is 5.3 x 10⁹ erythrocytes/ml of blood, the G6P-DH activity will result 695 / 5.3 = 131 mU/10⁹ erythrocytes.

To calculate the G6P-DH activity as mU/g hemoglobin, use the following formula:

$$G6P-DH \text{ (mU/g Hb)} = \frac{\text{mU erythrocytes per ml} \times 100}{\text{Hb (g/dl)}}$$

where:

100=factor to convert ml into dl

Hb (g/dl)=hemoglobin concentration for each sample

For example, if you obtain 695 mU for the erythrocytes present in 1 ml of blood, and the determined hemoglobin concentration is 15.0 g/dl, the G6P-DH activity will result (695 x 100) / 15.0 = 4633 mU/g Hb.

REFERENCE VALUES

In serum an activity under normal conditions is not virtually detectable (0-0.18 mU/ml).

In erythrocytes: 118 – 144 mU/10⁹ erythrocytes.

Each laboratory should define its own reference values.

PERFORMANCE CHARACTERISTICS

Linearity: up to 50 mU/ml in serum (at 37°C) and 1830 mU in erythrocytes for 1 ml of blood (25°C).

For higher values, dilute 1 volume of sample with 9 volumes of saline solution, repeat the assay and multiply the result by 10.

Within-run precision:

	Level 1	Level 2
Average (mU/10 ⁹ erythrocytes)	130	560
DS	1.10	8.82
CV %	0.85	1.57

Between-run precision:

	Level 1	Level 2
Average (mU/10 ⁹ erythrocytes)	108	502
DS	1.98	15.5
CV %	1.83	3.09

Correlation: G6P-DH FAR kit shows a correlation coefficient equal to 0.98, in comparison to another kit available on the market.

NOTES

- (*) dangerous reagent are marked with an asterisk. Refer to safety data sheet.
- Do not pipette by mouth. Apply normal precautions required to handle laboratory reagents.
- Chemistry analyzer parameters are available.

DISPOSAL

The product must be used for professional analysis only. The product must be disposed of according to national/international laws.

WARNINGS AND PRECAUTIONS



REAGENT 4: **H302** Harmful if swallowed.

REFERENCES

- A. Kornberg *et al.*, Methods in Enzymology I, Academic Press, New York, p. 323 (1955)

MANUFACTURER

FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

tel +39 045 6700870

website <http://www.farddiag.com>

e-mail: order@farddiag.com

e-mail: farddiag@farddiag.com

KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalog number
	temperature limits
	use by
	attenzione
	consult accompanying documents