

# NAG

Colorimetric Determination  
of N-Acetyl-β-Glucosaminidase  
in Urine and in Serum

4 x 10 Tests

REF CY08-40T

Available for quality control:

1 x 1 ml AAP NAG STANDARD

REF 7506

## PRINCIPLE

N-acetyl-β-glucosaminidase (NAG) catalyzes the hydrolysis of p-nitro-phenyl-N-acetyl-β-D-glucosaminide to p-nitrophenol and N-acetylgluco-samine. The liberated p-nitrophenol is proportional to the enzymatic activity and is colorimetric defined in an alkaline medium.

## REAGENTS

Kit components:	REF CY08-40T	Quantity
REAGENT 1 (lyo) Substrate	CY08-40TR1	4 vials
REAGENT 2 Buffer pH 5.2	CY08-40TR2	2x15 ml
(*) REAGENT 3 Alkalizer	CY08-40TR3	2x17 ml

STABILITY: stored at 2-8°C to keep reagents stable up to the expiration date on the label.

## PREPARATION OF REAGENT 1

REF 7114

Reconstitute the contents of a vial of Reagent 1 with exactly 5.2 ml of distilled water. Shake gently until complete dissolution.

STABILITY: 3 days at 2-8°C.

## SAMPLE

Non dialyzed urine.

STABILITY: 2 days at 2-8°C, 2 months at -20°C.

Serum: dilute the serum 1:4 with saline solution.

STABILITY: non diluted serum: 10 days at 2-8°C, 2 months at -20°C.

## MANUAL ASSAY PROCEDURE

Wavelength:	spectrophotometer 405 nm photometer 400-420 nm
Optical path:	1 cm
Reading:	against distilled water
Temperature:	37°C
Reaction time:	15 minutes
Linearity:	urine : up to 60 U/L serum: up to 240 U/L

Pipette into test tubes or cuvettes labeled as it follows:

S: sample, B/S: blank sample, B/R: blank reagent:

	S	B/S	B/R
Reagent 1	0.5 ml	---	0.5 ml
Reagent 2	---	0.5 ml	---

After reagents reach 37°C, add mixing accurately:

Sample	0.1 ml	0.1 ml	---
Saline solution	---	---	0.1 ml

Incubate for exactly 15 minutes at 37°C and add mixing accurately:

Reagent 3	0.400 ml	0.400 ml	0.400 ml
-----------	----------	----------	----------

Within 30 minutes, read the sample (As), the blank sample (Abs) and the blank reagent (Abr) absorbencies against distilled water.

## CALCULATION

URINE:

Activity in U/L = (As - Abs - Abr) x 36

Activity per mg of creatinine in mU/mg of creatinine =  
(U/L x 100) / mg of creatinine for 100 ml

SERUM:

Activity in U/L = (As - Abs - Abr) x 144

## REFERENCE VALUES

URINE:

Concentration: 5.7 ± 1.9 U/L

4.2 ± 1.2 mU/mg of creatinine

Excretion rate:

male 4.3 ± 1.3 mU/min

female 3.9 ± 1.3 mU/min

SERUM:

Concentration: 26 ± 3,1 U/L

## QUALITY CONTROL

All Clinical Chemistry laboratories should implement a quality control program. Control Urine of human origin are available for this purpose on request:

## PERFORMANCE CHARACTERISTICS

Linearity: up to 60 U/L in urine and to 240 U/L in serum.

For higher values, dilute 1 volume of sample with 9 volumes of saline solution, repeat the test and multiply the result by 10.

Within-run precision (urine sample):

Average (U/L)	Level 1	Level 2
DS	5.01	35.1
CV %	0.043	0.407
	0.86	1.16

Between-run precision (urine sample):

Average (U/L)	Level 1	Level 2
DS	5.50	32.5
CV %	0.085	0.828
	1.54	2.55

Correlation: FAR NAG kit shows a correlation coefficient of 0.98 in comparison to another kit available on the market.

## DISPOSAL

The product must be used for professional analysis only. The product must be disposed of according to national/international laws.

## WARNINGS AND PRECAUTIONS



Reagent 3 may cause severe eye irritation (H319). It also cause skin irritation (H315). In case of skin contact: rinse thoroughly with water. In case of contact with the eyes: rinse thoroughly for several minutes. If eye irritation continues, call a doctor.

## NOTES

- (\*) dangerous reagents are marked by an asterisk. Refer to safety data sheet.
- Centrifuge turbid samples at 3000 rpm x g for about 5 minutes.
- One blank reagent is enough for each assay series; if obtained values are constant, the determination can be done periodically.
- Urine pH and protease presence have a great influence on urinary enzyme stability. NAG is more stable with pH between 5 and 7. It is recommended to adjust the pH of highly acid or alkaline urines.
- Each laboratory should define its own reference values.
- Reaction volumes can be proportionally changed.
- Disposal waste according to local laws.
- Chemistry analyzer parameters are available.

## REFERENCES

- Maruhn D., Clin. Chim. Acta 73 (1976) 453
- Gressner A. M., Roebruck P., Clin. Chim. Acta 124 (1982) 315

## MANUFACTURER

FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

tel +39 045 6700870

website <http://www.fardiag.com>

e-mail: [order@fardiag.com](mailto:order@fardiag.com)

e-mail: [fardiag@fardiag.com](mailto:fardiag@fardiag.com)

## KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalog number
	temperature limits
	use by
	caution
	consult accompanying documents

Edition 01 - Jan 2021 RR