

# GLYCOSYLATED HEMOGLOBIN A<sub>1c</sub> – PRONTO

Chromatographic – spectrophotometric  
determination of Hemoglobin A<sub>1c</sub> in blood  
Ion exchange - Independent temperature method

20 tests

**REF** KR05-20

## INTENDED USE

Kit for quantitative *in vitro* determination of Hemoglobin A<sub>1c</sub> in blood.

## PRINCIPLE

By mixing a quote of whole blood with a hemolyzing reagent to obtain the lysis of the red blood cells is obtained, with also the release of hemoglobin and the elimination of the labile fraction.

Then hemoglobin A<sub>1c</sub> fraction is separated from HbA<sub>1A+B</sub> by cationic exchange chromatography.

The ratio between the eluated absorbance (at 415 nm) and the total hemoglobin absorbance provide the percentage of HbA<sub>1c</sub> in the sample. The procedure doesn't require any calibration and can be performed with a temperature range of 20-28°C.

## REAGENTS AND COLUMNS

Kit components:

**REF** KR05-20

**REAGENT 1** Potassium biphthalate 1 x 10 ml

**REAGENT 2** Buffer 42 mmol/L, LiCl, pH 5.96 1 x 40 ml

**REAGENT 3** Buffer 50 mmol/L, LiCl, pH 5.98 1 x 40 ml

**COLUMNS** Chromatographic columns 20

NOTE: use only reagents and columns of the same batch.

(\*) Dangerous reagents are marked by an asterisk. Refer to MSDS.

STABILITY: stored at 15-30°C, reagents and columns are stable up to the expiration date on the label. Store columns in the dark.

## SAMPLE

Whole blood collected with EDTA or oxalate fluoride. Do not use heparin as anticoagulant.

STABILITY: 7 days at 2-8°C.

## MANUAL ASSAY PROCEDURE

Wavelength: 415 nm  
Optical path: 1 cm  
Reading: against distilled water  
Temperature: 20-28°C  
Method: spectrophotometric  
C.V. (intra-assay): < 1.5%  
C.V. (inter-assay): < 4.0%

NOTE: Perform the assay with a temperature range between 20 and 28°C. Reagents and columns must have the same temperature before use.

## PREPARATION OF HEMOLYSED SOLUTION

Pipette in a tube:

Blood	0.050 ml
Reagent 1	0.450 ml

Mix vigorously and incubate at room temperature (20-28°C) for 10-15 minutes.

NOTE: use the hemolysed solution within 2 hours from its preparation.

## CHROMATOGRAPHIC SEPARATION

WARNING: during the whole chromatographic separation procedure, do not leave the resin without buffer for more than 5 minutes.

Pipette into the column:

Hemolysed solution	0.050 ml	Wait for 3 minutes. Discard any eventual eluate
--------------------	----------	---

Pipette into the column:

Reagent 2	2.0 ml	discard the eluate
-----------	--------	--------------------

Place the column in a tube (16 x 160 mm) and pipette:

Reagent 3	2.0 ml	collect the eluate (HbA <sub>1c</sub> fraction)
-----------	--------	---

Mix the obtained eluate and read the HbA<sub>1c</sub> fraction absorbance at 415 nm against distilled water (A Hb A<sub>1c</sub>).

## PREPARATION OF TOTAL HEMOGLOBIN

Pipette in a tube (16 x 160 mm):

Hemolyzed solution	0.100 ml
Distilled water	12.0 ml

Mix and read the absorbance of total hemoglobin (A Hb TOTAL) at 415 nm against distilled water.

## CALCULATION

$(A \text{ Hb A}_{1c}) / (3 \times A \text{ Hb TOTAL}) \times 100 = \% \text{ HbA}_{1c}$

## REFERENCE VALUES

Normal range: 4.2 - 6.2 %

These values are only indicative; each laboratory should define its own reference values.

## NOTES

- Interferences: wrong values may result from samples with abnormally high quantities of other hemoglobins due to their simultaneous elution with HbA<sub>1c</sub> (HbF) or to their differences in glycation compared with HbA's (HbS).
- The comparison between FAR kit (Y) with another kit (X) available on the market to determine A<sub>1C</sub> glycosilated hemoglobin, gave the following correlation line:  
 $Y=1.45X-3$   $R=0.957$
- Disposal of all waste material should be in accordance with local regulations.

## REFERENCE

- Mayer et Freedman (Clin.Chim.Acta, 1983; 127: 147-184)

## MANUFACTURER



FAR

Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY

phone +39 045 6700870

website <http://www.farddiag.com>

e-mail: [order@farddiag.com](mailto:order@farddiag.com)

e-mail: [farddiag@farddiag.com](mailto:farddiag@farddiag.com)

## KEY SYMBOLS

	In Vitro diagnostic medical device
	batch number
	catalogue number
	temperature limits
	use by
	caution
	read instructions for use

**IVD**



Ed. 02 - Jan 2024 RR