RID

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INTENDED USE

Kit for quantitative determination in radial immunodiffusion of human plasmaproteins in serum and plasma.

CLINICAL SIGNIFICANCE

Refer to the data reported in the bibliography.

PRINCIPLE

The sample (serum or plasma) containing the antigen (plasmaprotein), inoculated in the well of the plate, diffuses radially in the agarose gel, reacting with specific antibodies incorporated in agarose gel and forms a visible ring of immune complex. The diameter of the ring is proportional to the concentration of the protein in the sample.

SAMPLE

Plasma, fresh or stored at 2-8°C for fibrinogen e antithrombin III.
Serum, fresh or stored at 2-8°C for all other proteins.
Centrifuge turbid samples before use.
STABILITY: 2 days at 2-8°C, 30 days at -18°C.

REAGENTS

RID plates contain, in the layer of agarose gel, the specific antiserum to human plasmaprotein obtained by immunizing rabbits and goats.
STABILITY: RID plates, if stored at 2-8°C in their original pouch, can be used up to the expiration date on the label. Avoid heat sources and thermal shocks.
Store RID plates horizontally and upside down, to avoid water drop formation inside the well. Discard the plates in case of deterioration signs.
Unused wells can be used within one week from opening, keeping the plate at 2-8°C, sealed in its original pouch with a damp paper towel.
DO NOT FREEZE: do not store RID plates near the cooling coil of the refrigerator.

MATERIAL REQUIRED NOT SUPPLIED

Multiparametric calibrator.
Micropipette to accurately distribute 5 μ l of sample.
Magnifier lupe with a tenth of millimeter scale.

ASSAY PROCEDURE

Open the pouch, remove the lid and wait about 5-10 minutes at room temperature to allow any eventual condensation to evaporate.
Do not pipette any sample or calibrator if any humidity is present in the wells.
Pipette in the well 5 μ l of undiluted or diluted patient sample as indicate in the following tables:

Undiluted sample	
IgA	C4
IgG	C3c+C4
IgM	Transferrin
IgA+IgG+IgM	Antitrombin III
C3c	

Diluted sample 1:10 (1+9)	
Dilute with PBS or isotonic salty solution	
Ig/Light chain κ	Ig/Light chain κ +Ig/Light chain λ
Ig/Light chain λ	

Close the plate with the lid, incubate at room temperature (about 22-25°C) for about 20 minutes until the sample is completely absorbed into the gel.
Return the plate to the pouch, zip the pouch and incubate at room temperature in the upside-down pouch for 48 hours.
In order to prevent the gel from drying out, it is advisable to also put a damp paper towel in the pouch.

EVALUATION OF THE RESULTS

After 48 hours, keeping the plate upside down, measure the diameter of the sample precipitation ring, using a suitable measuring device (mm scale) with 0.1 mm accuracy.
In case of doubt regarding the complete diffusion, check the diameter after 24 hours, to make sure no diameter increasing happened.
Read on the conversion table the concentration value of the searched protein as expressed in mg/dL.

REFERENCE VALUES

See attached table of reference values.

LIMITS OF THE METHOD

The limits of the method are indicated on the conversion tables.
For rings with a diameter greater than the maximum indicated in the table, repeat the test with diluted sample and consider the dilution.

QUALITY CONTROL - CALIBRATION

Precision: CV < 4%

Accuracy: \pm 9% reference value.

Calibration: it is recommended to use calibrators to verify operating accuracy.
Deviations greater than \pm 0.2 mm from the value assigned to the calibrator are significant. They can depend on environmental conditions, error of the reading or measuring system and operator accuracy, for which they are systematic errors. In this case, take them into account when assigning the values to the samples.

DISPOSAL

After use, all material that came in contact with the sample must be considered potentially infected and disposed according to the law in force.

PRECAUTIONS

The plates contain sodium azide as preservative at a concentration <1%, therefore they are classified as not dangerous and do not require the safety data sheet. Do not ingest and avoid contact with skin.


BIBLIOGRAPHY

1. Parker C. W., Ed., Clinical Immunology, Saunders Company, Vol I,II (1980).
2. Putnam F. W., Ed. The plasma proteins, Academic Press, Vol 1,2 (1977).
3. Dammacco F., Ed., Ermes, Immunologia in medicina, (1989).

PRODUCER

FAR
Via Fermi, 12 - 37026 Pescantina - VERONA - ITALY
Phone +39 045 6700870
website: <http://www.farddiag.com>
e-mail: farddiag@farddiag.com / order@farddiag.com

SYMBOLS TABLE

	In Vitro diagnostic medical device
	batch number
	catalogue number
	temperature limits
	use by
	caution
	consult accompanying documents

Issue 01 - Jun. 2021 RR